

MicroActive® Resveratrol-SR

(Micronized Resveratrol in sustained release matrix)

Product Specifications

Raw Material: trans-Resveratrol 98%

Appearance: Off white granular powder

Ingredients: HPMC, Resveratrol, Carnauba wax, β -Cyclodextrin

Solubility: Soluble in water, methanol, ethanol

Quantitative Analysis: HPLC

Assay: Not less than 30% Resveratrol

Heavy Metals: <10ppm

Loss on Drying: Not more than 10%

Pesticide Residues (organochlorine, organophosphorus): None detected

Microbiological:

Total Plate Count: NMT 3000cfu/g

Salmonella sps, E.Coli, S. aureus, P. aeruginosa : Not detected

Yeast and Mold: NMT 100cfu/g

Storage: Store between 4-10°C protected from light, heat and moisture

Shelf-life: Two years from date of manufacture

Use: Suitable for capsules

